

Virucidal Cleaner Disinfectant Soluble Sachet

Cleanline® Eco Virucidal Cleaner Disinfectant soluble sachet is a multi-purpose cleaner, effective against bacteria, yeasts, enveloped viruses and conforms to EN1276, EN13697, and EN14476

- Fragrance-free, chlorine-free and alcohol-free cleaner disinfectant
- Proven effective against bacteria, yeasts, enveloped viruses
- Phosphate-free and contains biodegradable surfactants
- Independently tested to European Standard Test Method BS EN 1276, EN 13697 and EN14476
- Prepared solution has no health or environmental hazards in use
- No dosing equipment required and minimal storage needed
- Biodegradable and compostable packaging, no single-use plastic
- No chemical contact, eliminating risk of spillages and skin/eye contact
- Lightweight and easy to handle, reducing strain on the user
- Control costs and eliminate overdosing with an easy one-dose sachet
- Product has a low carbon footprint as water is added at the point of cleaning
- 95% reduced packaging compared to ready-to-use products
- Formulations are Vegan Society approved – no animal testing or animal-derived ingredients
- Minimal staff training required - The pre-dosed soluble sachet is easy to prepare and use, simply add 1 sachet to 600ml of warm water in a Cleanline **T11** Disinfectant reuse and refill bottle



CODE AND CASE SIZE:

CL4079 - 20 sachets

CL4079-S - Pack of 4 x 20 sachets **CL9011** - Empty

T11 labelled trigger spray bottle

PRODUCT INFORMATION:

UFI: 0RXA-G6YM-430P-5GJM

PCS No. TBC

Appearance: White

Odour: No odour

Oxidising: Non-oxidising

Solubility: Fully soluble in water

pH as supplied: 10-11 (@1%)

Packaging: Compostable and biodegradable recycled paper pouch

Contains: ALKYL (C12-14)

DIMETHYLBENZYLAMMONIUM CHLORIDE

ALCOHOLS C9-11, ETHOXYLATED

QUALITY ASSURANCE:

This product is manufactured in the UK for Bunzl UK & Ireland and produced under ISO 9001 Quality Management System & ISO 14001 Environmental Management Systems. This product has not been tested on animals.

ENVIRONMENTAL:

All surfactants used in Cleanline products comply with the current European regulations concerning biodegradability & protection of the environment.

SAFE HANDLING:

Always wear the correct Personal Protective Equipment (PPE) when handling, dispensing or using chemicals. All PPE should comply with appropriate standards and be readily available for use.



In case of an accident please refer to the product Safety Data Sheet (SDS) for immediate first aid measures and helpline telephone number for a medical emergency.

SAFE STORAGE:

Store the sachet in a dry area and keep them inside a sealed packet until needed

(UK) PrimeSource, PO BOX 15247,

Birmingham, B23 3HN, UK

(IE) PrimeSource, Unit D9, Horizon Logistics Park,

Swords, Co.Dublin, K67 N4T2, Ireland

Email: info@prime-source.co.uk

ISSUE DATE: 24/10/2022



Product information continued for Cleanline Eco Virucidal Cleaner Disinfectant Soluble Sachet

ENSURE YOU ARE WEARING RECOMMENDED PERSONAL PROTECTIVE EQUIPMENT (PPE) AND HAVE THE CORRECT EQUIPMENT FOR THIS TASK.

DIRECTIONS ON HOW TO PREPARE:



1. Add small amount of warm water to to T11 labelled spray bottle



2. Add 1 sachet, replace trigger head and gently shake bottle



3. Fill with warm water to 600ml mark



4. Replace trigger head and shake well

DIRECTIONS FOR USE:

Two stage cleaning (best practice)

Stage 1 Clean: Remove all loose dirt/soiling and food debris from surface. Clean the surface / equipment and agitate any stubborn soiling with a scouring pad to remove all visible soiling.

Stage 2 Disinfect: Apply the solution to the surface/equipment and allow recommended contact time (5 mins for protocol EN1276, EN14476 and EN13697) . Wipe with a clean cloth and rinse thoroughly with clean water. Allow to air dry.

IMPORTANT: IF ANY SACHETS ARE DAMAGED PLEASE DO NOT HANDLE THEM WITHOUT WEARING PPE. ONLY HANDLE THE SACHET WITH DRY HANDS.

MICROBIOLOGICAL ACTIVITY

Tested independently in accordance with;

BS EN 1276:2019

Tested under clean conditions, with a contact time of 5 minutes when product is diluted using 1 x 15g sachet per 750ml (1.3%v/v)

Bacterial action against *Pseudomonas aeruginosa* ATCC 15442, *Staphylococcus aureus* ATCC 6538, *Escherichia coli* ATCC 10536 and *Enterococcus hirae* ATCC 10541

BS EN 1650:2019

Tested under clean conditions, with a contact time of 5 minutes when product is diluted using 1 x 15g sachet per 750ml (1.3%v/v)

Microbial strains in the test: *Candida albicans* ATCC 10231

BS EN 13697:2015+A1:2019

Tested under clean conditions with a contact time of 5 minutes when product is diluted using 1 x 15g sachet per 750ml (1.3%v/v)

Bacterial action against *Pseudomonas aeruginosa* ATCC 15442, *Escherichia coli* ATCC 10536, *Enterococcus hirae* ATCC 10541, *Staphylococcus aureus* ATCC 6538

A summary of the tests is available on request.



CODE AND CASE SIZE:

CL4079 - 20 sachets

CL4079-S - Pack of 4 x 20 sachets

CL9011 - Empty **T11** labelled trigger spray bottle

PRODUCT INFORMATION:

UFI: 0RXA-G6YM-430P-5GJM

PCS No. TBC

Appearance: White

Odour: No odour

Oxidising: Non-oxidising

Solubility: Fully soluble in water

pH as supplied: 10-11 (@1%)

Packaging: Compostable and biodegradable recycled paper pouch

Contains: ALKYL (C12-14)

DIMETHYLBENZYLAMMONIUM CHLORIDE

ALCOHOLS C9-11, ETHOXYLATED

QUALITY ASSURANCE:

This product is manufactured in the UK for Bunzl UK & Ireland and produced under ISO 9001 Quality Management System & ISO 14001 Environmental Management Systems. This product has not been tested on animals.

ENVIRONMENTAL:

All surfactants used in Cleanline products comply with the current European regulations concerning biodegradability & protection of the environment.

SAFE HANDLING:

Always wear the correct Personal Protective Equipment (PPE) when handling, dispensing or using chemicals. All PPE should comply with appropriate standards and be readily available for use.



In case of an accident please refer to the product Safety Data Sheet (SDS) for immediate first aid measures and helpline telephone number for a medical emergency.

SAFE STORAGE:

Store the sachet in a dry area and keep them inside a sealed packet until needed

(UK) PrimeSource, PO BOX 15247,

Birmingham, B23 3HN, UK

(IE) PrimeSource, Unit D9, Horizon Logistics Park,

Swords, Co.Dublin, K67 N4T2, Ireland

Email: info@prime-source.co.uk

ISSUE DATE: 24/10/2022



Cleanline® is an exclusive brand of Bunzl UK & Ireland

© 2022 Bunzl UK Ltd. All Rights Reserved.